




Fever thermometers

CE Certified Medical Devices
(E.D. 93/42 / EEC)



- Fevertest**
Certified Medical Device
- speed of response (15 seconds)
 - reusable
 - pocket size
 - unbreakable
 - no batteries needed

Fevertest - Certified Medical Device




- Fevercontrol**
Certified Medical Device
- fever test in 15 seconds
 - easy reading
 - reusable
 - long lasting

Fevercontrol - Certified Medical Device



- 48h Fevermonitor**
Certified Medical Device
- designed for babies and children
 - monitor the fever status for up to 48 hours
 - even during night
 - test in 30 seconds
 - unbreakable, reliable and comfortable



48h Fevermonitor - Certified Medical Device



International Products & Services

Baby care

Thermometers for babies




bathermo
termometro da bagno - bath thermometer
thermomètre de bain - termómetro de baño

infrangibile
lunga durata
preciso

unbreakable
long lasting
precise

incassabile
durabile
preciso

irrompibile
lunga duración
exacto

Thermometer which allows to control the baby bath water in a few seconds.

Bathermo




idealtemp
termometro per biberon
feeding bottle thermometer
thermomètre pour le biberon

universale - riutilizzabile - lunga durata
universal - reusable - long lasting
universel - réutilisable - longue durée

Unbreakable thermometer that shows the ideal feeding temperature.

Idealtemp



Beauty Kit

Cosmetic testers for personal use

100% made in italy

BEAUTY KIT is a unique and innovative product for facial skin care. Beauty kit includes 3 practical disposable cosmetic testers:

- SKIN TEST: skin type indicator.
- HYDRATION TEST: skin hydration level indicator.
- WRINKLE TEST: facial expression lines and wrinkle indicator.

The use is simple and the result is fast. With BEAUTY KIT it is possible to carry out tests independently for a complete skin check-up and to choose the most suitable cosmetic treatment.

• Disposable
• Non-toxic
• Harmless

BeautyKit



Perfumed patches

With natural microencapsulated essential oils

- Made with natural essential oils unpleasant to mosquitoes and moths
- Non-toxic
- Suitable for babies and children
- Long lasting
- Reactivable



1. Peel the patch from the box.
2. Stick the patch to the surface.
3. Peel the patch from the surface.
4. Reuse the patch.

Perfumed patches

